

DYNAMIC Duo

Check out these impressive young sires recently added to the Accelerated Genetics Jersey lineup!



JX Rowleys 1239 Leonel 1749 {4} (E-90%), grandam of JX MAYPOP {5}-P



JX Victory 50891 Chief 56815 {5} (VG-85%), grandam of HOVLAND

614JE2235

JX MAYPOP {5}-P

MONSTER {6}-P x JX STONEY {3} x JX LEONEL {3}
BBR 100 A2A2 K-Casein: BB JH1F JNSC

+179	+16.0	+115	+0.08%	+0.15%
GJPI™	JUI™	CFP	Protein	Fat
+787	+993	2.83	105	
HHP\$®	DWP\$®	SCS	Z Mast	



JX MAYPOP {5}-P combines exceptional mastitis resistance and solid components with strong calf health traits.

614JE2240

HOVLAND

UNCOMMON x WATKINS x JX Chief {6}
BBR 100 A2A2 K-Casein: BB JH1F JNS-TF

+176	+14.4	+1,052	+114	+707
GJPI	JUI	Milk	CFP	CM\$
+778	+926	+4.6	+106	
HHP\$	DWP\$	PL	Z Mast	



HOVLAND is your source for elite production, valuable Type, and fitness trait improvement.

Contact your local representative to add this impressive duo to your herd today!

614-873-4683 | WWW.ACCELGEN.COM



02/24 Zoetis data. 12/23 CDCB-SAJICA Genomic Evaluation Ref. %: JX MAYPOP {5}-P Type 78, Yield 75, SCS 74, Z Mast 35; HOVLAND Type 78, Yield 74, PL 71, Z Mast 35. *JPI and JUI are trademarks of the American Jersey Cattle Association. gender SELECTED is a trademark of Select Sires Inc. ®DWP\$ is a registered trademark of Zoetis Inc., its affiliates and/or its licensors. Mastitis ResistantPRO, FeedPRO, RobotPRO and HHP\$ are registered trademarks of Select Sires Inc. JX MAYPOP {5}-P and HOVLAND both qualify for semen export to Canada and are only available sexed. All gender SELECTED semen is processed using Ultraplus™ technology. Ultraplus is a trademark of STGen LLC. Limitations and Conditions of Sale: gender SELECTED semen shall be used only for the single insemination of one female bovine during natural ovulation with the intent to produce a single offspring unless specifically approved on an individual customer basis by Sexing Technologies in writing. As a condition of purchasing gender SELECTED semen, the purchaser agrees that gender SELECTED semen will not be reverse sorted or re-sorted unless specifically permitted, in advance, on a case-by-case basis by Select Sires in writing. Select Sires intends to monitor the use of the gender SELECTED semen and vigorously enforce these restrictions on use. Please see <http://www.selectsires.com/designations/gendersselected.html> for additional details.

= genderSELECTED* = Mastitis ResistantPRO* = FeedPRO* = RobotPRO*